

## Introduction

The DogCatcher Nanoselector® Magnetic beads consists of DogCatcher protein, which is covalently bound to magnetic beads. DogCatcher Nanoselector® Magnetic beads is used to immunoprecipitate DogTag-fusion proteins from cell extracts of various organisms like mammals, plants, bacteria, yeast, insects etc.

## Properties

<b>Ligand</b>	DogCatcher
<b>Reactivity</b>	Specifically binds to DogTag and its fusion protein
<b>Binding capacity</b>	20 µg of recombinant DogTag-GFP per 20 µL bead slurry
<b>Bead size</b>	~ 40 µm
<b>Storage buffer</b>	10mM PBS (pH 7.5), 0.01% proclin300
<b>Storage conditions</b>	Upon receipt store at +4°C. Do not freeze!
<b>Stability</b>	Stable for 1 year upon receipt.
<b>Shipment</b>	Shipped at low temperature(2-8°C)

## Suggested buffer compositions

### Required buffer solutions

Buffer	Composition
Lysis buffer	10 mM Tris-HCL pH 7.5, 150 mM NaCl, 0.5 mM EDTA, 0.5 % NP40
RIPA buffer	10 mM Tris-HCL pH 7.5, 150 mM NaCl, 0.5 mM EDTA, 0.1 % SDS, 1 % TritonX-100, 1 % deoxycholate
Dilution buffer	10 mM Tris-HCL pH 7.5, 150 mM NaCl, 0.5 mM EDTA
Wash buffer	10 mM Tris-HCL pH 7.5, 150 mM NaCl, 0.05 % NP40, 0.5mM EDTA
2xSDS Sample buffer	120 mM Tris-HCL pH 6.8, 20 % glycerol, 4 % SDS, 0.04 % bromophenol blue, 10 % β-mercaptoethanol
Acidic elution buffer	200 mM glycine pH 2.5
Neutralization buffer	1 M Tris-HCL pH 10.4

Note: Use your equivalent cell lysis buffer for other cell types like yeast, plants, insects, bacteria.

Note: Consider using a Wash buffer without detergent for co-immunoprecipitation.

# DogCatcher Magnetic beads

Code: 086-103-003



## Product sizes

Product	Size
DogCatcher Nanoselector® Magnetic beads	20 reactions (500 µL slurry)
	200 reactions (5 mL slurry)

## Immunoprecipitation protocol

### Cell material

The following protocol describes the preparation of mammalian cell lysate! For other type of cells, we recommend using 500 µg of cell extract and start the protocol with step Bead equilibration.

### Mammalian cell lysis

Note: Harvesting of cells and cell lysis should be performed with ice-cold buffers. We strongly recommend to add protease inhibitors to the Lysis buffer to prevent degradation of your target protein and its binding partners. For one immunoprecipitation reaction, we recommend using  $\sim 10^6$  -  $10^7$  cell.

1. Choice of lysis buffer: · For cytoplasmic proteins, resuspend the cell pellet in 200 µL ice-cold Lysis buffer by pipetting up and down. Supplement Lysis buffer with protease inhibitor cocktail and 1 mM PMSF (not included). · For nuclear/chromatin proteins, resuspend cell pellet in 200 µL ice-cold RIPA buffer supplemented with DNaseI (f.c. 75-150 Kunitz U/mL), MgCl<sub>2</sub> (f.c. 2.5 mM), protease inhibitor cocktail and PMSF (f.c. 1 mM) (not included).
2. Place the tube on ice for 30 min and extensively pipette the suspension every 10 min.
3. Centrifuge cell lysate at 17,000x g for 10 min at +4°C. Transfer cleared lysate (supernatant) to a precooled tube and add 300 µL Dilution buffer supplemented with 1 mM PMSF and protease inhibitor cocktail (not included). If required, save 50 µL of diluted lysate for further analysis (input fraction).

### Plant tissue lysis (e.g. Arabidopsis thaliana):

1. Use a tissue to absorb the liquid cultured Arabidopsis seedlings, weigh 0.7 grams, and transfer them to a mortar.
2. Freeze with liquid nitrogen and grind into fine powder with a pestle.
3. Add 4ml of extraction buffer (PBS, 0.5% Triton X-100, 0.5 mM PMSF, protease inhibitor), slowly thaw, and further grind the sample.
4. Transfer the sample to a 1.5 ml microcentrifuge tube and centrifuge at 16000 x g for 20 minutes at 4 ° C.
5. Pass the supernatant through a 0.20 µ m filter syringe and place the extract on ice. If required, save 50 µL of diluted lysate for further analysis (input fraction).

## Bead equilibration

1. Resuspend the beads by gently pipetting up and down or by inverting the tube. Do not vortex the beads!
2. Transfer 25  $\mu$ L of bead slurry into a 1.5 mL reaction tube.
3. Add 500  $\mu$ L ice-cold Dilution buffer.
4. Separate the beads with a magnet until the supernatant is clear. Discard the supernatant.

## Protein binding

1. Add diluted lysate to the equilibrated beads.
2. Rotate end-over-end for 1 hour at +4°C.

## Washing

1. Separate the beads with a magnet until the supernatant is clear.
2. If required, save 50  $\mu$ L of supernatant for further analysis (flow-through/non-bound fraction).
3. Discard remaining supernatant.
4. Resuspend beads in 500  $\mu$ L Wash buffer.
5. Separate the beads with a magnet until the supernatant is clear. Discard the remaining supernatant.
6. Repeat this step at least twice.
7. During the last washing step, transfer the beads to a new tube. Optional: To increase stringency of the Wash buffer, test various salt concentrations e.g. 150-500 mM, and/or add a non-ionic detergent e.g. Triton X-100.

## Nanoselector related products

Code Number	Product Description	Applications	Size
002-101-002	Myc tag Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
002-101-003	Myc tag Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
019-101-002	GFP Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
019-101-003	GFP Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
016-101-002	DYKDDDDK tag Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
016-101-003	DYKDDDDK tag Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
003-101-002	HA tag Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
003-101-003	HA tag Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
015-101-002	MBP Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
015-101-003	MBP Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
010-101-002	GST Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
010-101-003	GST Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
011-101-002	SNAP tag Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
012-101-002	Halo Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
013-101-002	mNeongreen Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
014-101-002	TurboGFP Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
004-101-002	c-His tag Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
004-101-003	c-His tag Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
049-101-002	mWasabi Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
017-101-002	TagFP Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
025-101-002	Rabbit IgG Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
025-101-003	Rabbit IgG Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
001-101-002	Mouse IgG Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
001-101-003	Mouse IgG Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
001-300-002	Mouse IgG&Rabbit IgG Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
001-300-003	Mouse IgG&Rabbit IgG Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
067-101-003	Streptavidin Magnetic beads	IP,CHIP,MS,Purification	1mL(40 rxns)
064-101-002	V5 tag Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
064-101-003	V5 tag Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)

## Nanoselector related products

Code Number	Product Description	Applications	Size
086-101-003	Spycatcher Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
086-101-002	Spycatcher Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
082-101-002	StayGold/mBaojin Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
082-101-003	StayGold/mBaojin Nanoselector Magnetic	IP,CHIP,MS,Purification	0.5mL(20 rxns)
020-101-002	RFP Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
020-101-003	RFP Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
023-101-003	Human IgG Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
026-101-002	tdTomato Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
026-101-003	tdTomato Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
028-101-002	mScarlet Nanoselector Agarose	IP,CHIP,MS,Purification	0.5mL(20 rxns)
028-101-003	mScarlet Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
078-101-003	Biotin Nanoselector Magnetic beads	IP,CHIP,MS,Purification	0.5mL(20 rxns)
100-100-100	Binding Control Nanoselector Agarose	IP,CHIP,MS,Purification	1mL(40 rxns)
100-100-200	Binding Control Magnetic beads	IP,CHIP,MS,Purification	1mL(40 rxns)